

# detect-a-chip from The Pet Chip Co Ltd



5 Hereward Road  
Eastbourne.  
BN23 6TG.  
Tel/Fax: 01323 469106  
Mobile: 07876 282233  
email: [zoochip@tlc1.co.uk](mailto:zoochip@tlc1.co.uk)  
Website: [www.zoochip.com](http://www.zoochip.com)

## detect-a-chip .... it reads

**all ISO Transponders, Trovan's ID-100 'Unique' & Avid encrypted**

**detect-a-chip** has been specifically designed to meet the needs of those who require an easy-to-use, 'pocket sized' microchip reader on a limited budget. This scanner's exceptional reading ability makes **detect-a-chip** ideal for 'use in the field', in Kennels and catteries, in zoos, animal quarantines, grooming parlours, vets, animal sanctuaries and Rescue organisations.

At time of going to print:

**detect-a-chip** reads all known *animal implantable microchips including:*  
*ALL FDX-A & FDX-B, ISO 11784/5 compliant transponders;*  
*FECAVA transponders; Trovan ID-100 'Unique'; Avid Encrypted transponders.*



## Performance data:

Dimensions: 5 inches x 2.66 inches x 1 inch.L x W x H.  
(metric: 125mm x 68mm x 25mm)  
Weight: 100gms  
Power supply: 9V Battery  
Read distance: 2 - 4 inches (metric 5 - 10 cms)  
depending on transponder type.  
One button menu  
Belt clip  
Automatically turns off when not in use.

## Price

Protect your scanner against those occasional 'drop' accidents with

## protect-a-case

This strong, durable and inexpensive rubberised sleeve *could* save users the cost of a new scanner.

**protect-a-case** comes in Red or Dark Blue.



## Price